



arco-matic® – The gynaecological examination and treatment chair

An award-winning design and state-of-the-art technology combined at excellent value for money.

The arco-matic® examination and treatment chair, with a width of upholstery of up to 680 mm at the back section and an individually adjustable head cushion with magnetic attachment, guarantees patients a pleasant, comfortable and safe sitting experience.

With an access height of only 550 mm, even older or physically challenged patients can easily take a seat on the chair, and even at its lowest level, the seat position remains horizontal.

Designed for a working load of 250 kg, the chair offers great safety and many treatment options even for adipose patients.

The ambient lighting and the SCHMITZ colour concept combine to create a pleasant and anxiety-reducing atmosphere.

SCHMITZ – Your partner for gynaecology, video colposcopy, urology and proctology.



reddot award
winner

Made in Germany



arco-matic®
product configurator

Configure your new arco-matic® now:
www.schmitz-medical.com

Key Features



Flat position, e.g. for ultrasound scans, or shock position in an emergency



Work flexibly and save time with four freely selectable memory positions.



The deep gyn. cut-out offers improved access to the examination area.



Flexible use and effective cleaning – even under the arco-matic®! (castors Ø 75 mm or Ø 100 mm)



Side rails at the seat section for attaching accessories, e.g. an LED examination light



The LED examination light provides a homogeneous light field with 30,000 Lux.



Pelvic elevation and electrically adjustable leg positioning systems ensure the ideal lithotomy position.



Matched to your treatment room, the ambient LED illumination creates a pleasant atmosphere.



High safe working load of 250 kg



The removable rinsing basin with convenient holder always remains in a horizontal position.



Compatible with a wide variety of leg positioning systems for attachment to the foot support mounts or to side rails




Side rails at the back section allow accessories such as arm rests and an infusion pole to be used.




arco-matic® combines functionality and design to form the central element of any modern consulting room.


Highest position



Overview




Lowest position



Maximum back section tilt 52°

Shock position



Trendelenburg 29°

Note: The heights given above do not take into account a mobile version.

