

Suspension arm system



Why choose MZ Liberec?



Design - Maintenance - Flexibility


- Simple cleaning and disinfection thanks to the uniform surface
- Antibacterial surface treatment technology
- Compatible with a broad range of accessories



Functions and features

Ceiling pendants with compact, spring-balanced arms serve as a broad range of applications in the healthcare industry and in other service industries.

They can support operating lights, cameras, monitors, security X-ray screens and many other specialized devices.



Technical features



Thanks to the rotating contacts, the arms can easily rotated 360°




Spring arms with easy adjustment of balance, range and stop function



Customer solutions - variability

SPR offers assorted mounting options



ceiling-mounted




wall-mounted

Options - SPR


Ceiling versions - SPR 10

SPR 10-1




Loading capacity of 1 arm - Max. 28 kg
Weight - 30 kg

SPR 10-2
- low ceiling solution



Loading capacity of 1 arm - Max. 28 kg
Weight - 42 kg

SPR 10-3



Loading capacity of 1 arm - Max. 28 kg
Weight - 55 kg

SPR 10-4




Loading capacity of 1 arm - Max. 28 kg
Weight - 65 kg

Options - SPR


Ceiling versions - SPR 11

SPR 11-1



Loading capacity of 1 arm - Max. 12 kg
Weight - 18 kg


SPR 11-2



Loading capacity of 1 arm - Max. 12 kg
Weight - 28 kg

SPR 23

Innovation



Loading capacity of 1 arm - Max. 4 kg
Weight - 15 kg


Options - SPR

Wall version


SPR 23-1W




SPR 10-1W



SPR 11-1W



SPR 10 with electrical terminal



Mobile version

SPZ 13

The mobile stand is rated for medical devices with a maximum weight of 5 kg.

The device is primarily used to hold a heat source used in neonatal departments in healthcare facilities.



Options SPR	SPR 10	SPR 11	SPR 23	SPZ 13
Length of rotating arm	500 - 1 100 mm	700 - 1 100 mm	700 - 1 100 mm	
Length of tilting arm	960 mm	800 mm	550 - 860 mm	480 mm
Maximum load	28 kg	12 kg	4 kg	5 kg
Ceiling height options	normal / low ceiling			
Colour options	RAL 9002, RAL 9010, RAL 9016, OEM – customer solution			
Design materiálu	Powder-coated aluminium construction with antibacterial treatment			



