

Superior Technology Incredible Performance **F15 Series**

Fetal & Maternal monitor




Evolutionary Generation of Fetal Monitoring

As an ultimate CTG machine, F15 series fetal and maternal monitor is the perfect combination of advanced technology and sophisticated algorithm for Obstetrics' new era. With state-of-art design, elegant appearance and strong mechanical structure, F15 Series is dedicated to bring the most efficient and innovative work flow to clinics.




Sophisticated CTG Alarm System

The essential purpose of fetal monitoring is to recognize abnormal conditions in time. Integrated with international CTG interpretation standard, F15 will issue an alarm when dangerous FHR waveform is recognized.



Sinusoidal FHR



Low Variability FHR

Innovative Wireless Application

Offering comprehensive wireless solutions to clinics, freeing midwives and mom from messy wires. Additionally, F15 Air supports non-stop monitoring when transferring mom-to-be, data can be transferred from one device to another seamlessly.



FT20

Integrated Maternal Parameters Module

- Support NIBP, SpO₂ and PR
- Wearable device provides mom with more choices with their ambulation and position




FECG-TA

As an advanced module for F15 Air, FECG-TA can acquire FHR and MHR from advanced surface electrodes system using FECG and MECG technology, and obtain EHG simultaneously.

With all in one design, moms and caregivers are bestowed with more choices when monitoring.

Accurate


- Autonomous MECG channel for calibration
- Precise R-R interval identification from ECG
- Reliable for women of high BMI



Extraction

All in one design

- FECG, MECG, FHR, MHR, UA and MFM can be obtained with one module.



Calculation



Comfort

- Cover all fetal positions
- Weight less than 150g
- Passive acquisition and non-invasive



Separation

F15 Air

Comprehensive wireless monitoring solution

All fetal and maternal parameters can be monitored wirelessly, empowering mom and clinician with more choices, providing a more comfortable experience throughout the labor process.





F15 F15 Air