

• Enteral Feeding System




Enteral giving set (spike)



Enteral giving set (bag)



Enteral syringe



Enteral feeding catheter



Easy to feed Safe to use
Enteral Feeding Pump
EP-60 Series



MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.



Add.: 12th Floor, Baiwang Research Building, 5158 West Shahe Road, Xili, Nanshan, 518055 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA
Tel: +86 755 26953369 E-mail: info@medcaptain.com
Website: www.medcaptain.com
MEDCAPTAIN is a trade mark of Medcaptain Medical Technology Co., Ltd.
Specifications are subject to change without prior notice.
©2023 Medcaptain Medical Technology Co., Ltd. All rights reserved.
EN-EP-60-V2.8

Distributor:




Enteral Feeding System

Integrated enteral feeding system
We provide enteral giving sets,
enteral feeding catheters and
enteral syringes



Auto Anti Occlusion

Even between two meals,
feeds will not clot



4.0" Touch Screen

Touch to set parameters
Quick training and easy operation
Two user interface themes
Teaching videos inside




Robust Housing

High ingress protection level IP34
Housing protects pump
in case of feed solution leaks



Optional Heater

Heater works and stops simultaneously with the pump
Thermal insulation and safe voltage
Without additional power cord
Can reduce irritations in gastrointestinal tract



Compact Design:

Weight: 1.2 kg
Unique clamp allows
install and removal within 1 second
Anti-skid hidden handle

Superb Parameters:

Accuracy: $\leq \pm 5\%$ Delivery rate: 1-1200 ml/h
USB upload and download

Strong Battery:

Charge for 4 hours, work for 22 hours

Smart Prime:

Fast prime at 2000 ml/h
Automatically stops when prime is finished

Exclusive Design Because we care



Heater

Hidden Handle

Clamp