



moments to care



Ambulatory Infusion System



MEDCAPTAIN MEDICAL TECHNOLOGY CO.,LTD.



Room 701, Block A, Building7, Shenzhen International Innovation Valley,
Dashi First Road, Nanshan, 518055 Shenzhen, Guangdong,
PEOPLE'S REPUBLIC OF CHINA

Tel: +86 755 26953369 E-mail: info@medcaptain.com

Website: www.medcaptain.com

MEDCAPTAIN is a trade mark of Medcaptain Medical Technology Co., Ltd.
Specifications are subject to change without prior notice.

©2024 Medcaptain Medical Technology Co., Ltd. All rights reserved.

EN-Ambulatory Infusion System -8P-Version 2.0


DISTRIBUTOR:



ELFensor P series Ambulatory Infusion Pump


Designed for precise and quantitative control of medications administered to the patient. Comply IEC60601-1-11 for home healthcare environment.

Professional for more scenarios




PIEB mode including (Programmed Intermittent Epidural Bolus)

- PIEB dose
- PIEB dose + PCEA dose
- PIEB dose + continuous dose



PCA mode including (Patient Controlled Analgesia)

- Continuous dose
- PCA dose
- Clinician bolus



Trapezia mode

Used when the treatment involves a slow and gradual increase in doses at the beginning and a gradual decrease in doses at the end. For example, total parenteral nutrition (TPN) support.



Connectivity

iCMS supports connectivity to any hospital information system or clinical information system via HL7. iCMS monitored all pumps treatment status.




CQI Report

CQI (continuous quality improvement) helps us improve infusion device management continuously.

DoseControl

DoseControl module allows healthcare professional teams to design drug library for ambulatory pumps.





Medication Cassette Reservoir

The medication cassette reservoir consists of a cassette assembly and an extension set assembly.




Easy to add medicine



Easy to Install



Ease of observation




Colour-coded options




(Lightproof optional)



Material optional		
PVC(TOTM)	TPU	TPU+PVC(TOTM)



Extension Set with check valve



Extension Set with Luer or NRFit Connector



Maximum duration of use:
7 days



Latex free




DEHP free

Administration Sets

Anti-free-flow protection

The Administration Sets are designed to be used with Medcaptain ELfensor pump to allow fluid delivery from an IV bag. The administration sets vary in length, configuration, and storage volume.



Administration Sets With Vented Spike

Model	Length (cm)	PV (ml)	Description	Pcs/Carton
GEN01002	300	8.2	Vented Spike, Clear, Male Luer Lock Connector With Check Valve, 0.2µm In-Line Filter	50
GEN02012	300	7.2	Vented Spike, Clear, Male Luer Lock Connector With Check Valve, 1.2µm In-Line Filter	50
PCA01002	300	8.1	Vented Spike, Blue Striped, Male Luer Lock Connector With Check Valve, 0.2µm In-Line Filter	50
EPI01002	300	7.6	Vented Spike, Yellow Striped, Male NRFIT Connector , 0.2µm In-Line Filter	50
EPI02002	300	8.0	Vented Spike, Yellow Striped, Rotating Male Luer Lock , 0.2µm In-Line Filter, In-Line check valve	50
PCA03002	300	6.8	Vented Spike, Light Resistant, Male Luer Lock Connector With Check Valve, 0.2µm In-Line Filter	50
EPI03001	300	7.8	Vented Spike, Yellow Striped, Male NRFIT Connector, 0.2µm In-Line Filter, In-Line check valve	50
PCA04002	300	7.4	Vented Spike, Clear, Male Luer Lock Connector With Check Valve, 0.2µm In-Line Filter	50

Administration Sets With Vented Spike and Y-port

Model	Length (cm)	PV (ml)	Description	Pcs/Carton
GEN03012	300	7.4	Vented Spike, Clear, Male Luer Lock Connector With Check Valve, 1.2µm In-Line Filter, Y-port(Needle free)	50

PCA02002	300	8.6	Vented Spike, Blue Striped, Rotating male luer lock, 0.2µm In-Line Filter, Y-port(Female luer lock connector)	50
TPN01012	300	7.2	Vented spike, clear, male luer lock connector with check valve, 1.2µm In-Line Filter, Y-port(Needle free)	50

Administration Sets With Drip Chamber Vented

Model	Length (cm)	PV (ml)	Description	Pcs/Carton
GEN110000	240	9.0	Drip chamber vented, Clear, Male Luer Lock Connector With Check Valve	50
GEN120000	240	12.2	Drip chamber vented, Clear, Male Luer Lock Connector With Check Valve, Y-port(Needle free)	50
GEN 130000	240	9.2	Drip chamber vented, Light Resistant, Male Luer Lock Connector With Check Valve, Y-port(Needle free)	50
ONC110000	240	9.2	Drip chamber vented, Light Resistant, Rotating Male Luer Lock , Y-port(Needle free)	50
ONC120000	240	10.8	Drip chamber vented, Clear, Rotating Male Luer Lock , Y-port(Needle free), In-line check valve	50
ONC130020	240	12.6	Drip chamber vented, Clear, Rotating Male Luer Lock , 0.2µm In-Line Filter, Y-port(Needle free) , In-line check valve	50

Administration Sets With Female luer lock

Model	Length (cm)	PV (ml)	Description	Pcs/Carton
GEN21002	240	7.2	Female luer lock Connector, Clear, Male Luer Lock Connector With Check Valve, 0.2µm In-Line Filter	50
PCA21002	240	6.6	Female Luer Lock Connector, Blue Striped, Male Luer Lock Connector With Check Valve, 0.2µm In-Line Filter	50
PCA22002	240	7.0	Female luer Lock connector, Blue Striped, Rotating Male Luer Lock, 0.2µm In-line Filter, Y-port (Female luer lock connector) , In-line check valve	50
PCA31002	240	7.0	Male Luer Lock Connector, Clear, Male Luer Lock Connector With Check Valve, 0.2µm In-Line Filter	50